COVID-19 (SARS-CoV-2) Antigen Test Midstream - Saliva

Single Use Only.
For in vitro diagnostic use only.
Please read the instruction carefully before use.


Operational Use Video


[Intended use]

This product is used for in vitro qualitative detection of the SARS-CoV-2 antigen in human saliva specimen, which is used to assist in the diagnosis of COVID-19. This product is intended for home self-testing as a rapid test for novel coronavirus infection. Both people who have close contact with COVID-19 infected patients and symptomatic people can be tested. However, preliminary results need to be confirmed based on clinical diagnostic results. Please do not make a medical relevance decision without consulting with your doctor. Users aged 10 and over can test by themselves. Users under 10 years old should be tested with supervision or assistance of an adult. Test within the first 7 days of symptom onset when viral shedding/viral load is at its highest.

[Materials and Components]


[Preparation before the test]


[Test Procedure]

NOTE: Please keep the temperature at 15 ~ 30 ℃ and the humidity at 20%-80% during the whole test.

[Interpretation of test results]

[Summary]

The novel coronavirus belongs to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection, asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main Manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases. Once infected with the SARS-CoV-2 virus, you may be hospitalized and some complications may occur. If without prompt treatment it may even lead to death.

[Test principle]

This product uses the double antibody sandwich method to detect the SARS-CoV-2 N protein. When the sample contains the coronavirus antigen, both the test line (T) and the control line (C) will appear, and the result will be positive. When the sample does not contain the coronavirus antigen or no coronavirus antigen is detected, the test line (T) will not appear, only control line (C) will appear.

[Limitations of inspection methods]

- This test kit is only used for in vitro diagnosis.
- This test kit is only used to detect COVID-19 antigen in human saliva samples. The results of other specimens may be wrong.

- This test kit is only used for qualitative detection and cannot indicate the level of SARS-CoV-2 antigen in the specimen.
- This test kit is only a clinical auxiliary diagnostic tool. If the result is positive, it is recommended to use other methods for further examination in time and the doctor's diagnosis shall prevail.
- This test does not determine the aetiology of the respiratory infection caused by micro-organisms other than the SARS-CoV-2 virus.
- This test can detect both the viable and the non-viable SARS-CoV-2 virus. False negative results may be given following poor sampling.
- Any failure to respect the test procedure may negatively impact the performance of the test and/or invalidate the test result.
- If the result of the test is negative, yet clinical symptoms persist, it is advised that you carry out additional tests using other clinical methods. A negative result at no time rules out the presence of antigens of the SARS-CoV-2 virus in the sample, as they may be present but at a level inferior to the minimum detection level of the test, or if the sample has been collected incorrectly.
- A negative result does not rule out infection by the SARS-CoV-2 virus, particularly in people who have come into contact with the virus. Follow-up tests with molecular diagnostics should be scheduled to rule out infection in these people. Persons who show symptoms of the disease but have a negative result until infection is ruled out should follow country-specific restrictions.
- This test is not a substitute for a medical consultation, or for the result of a biological analysis carried out in a medical analysis laboratory.
- Positive test results do not exclude the possibility of co-infections of other pathogens.

[Warnings and Precautions]

- Read the instructions carefully before using the kit, and strictly control the reaction time. If you do not follow the instructions, you may get inaccurate results.
- Do not eat, drink, chew gum, smoke or vape for at least 30 minutes before testing. False negative results may occur if the saliva is not collected properly.
- Guard against moisture, do not open the foil pouch before it is ready for testing. Do not use it if the foil pouch is damaged or the test device is damp.
- Please use it within the validity period.
- Do not replace the components in this kit with components in other kits.
- The kit shall be stored in strict accordance with the conditions specified in this manual. Please do not store the kit under freezing conditions.
- The test methods and results must be interpreted in strict accordance with this specification.
- Negative results may occur if the SARS-CoV-2 antigen titer in the specimen falls below the minimum detection limit of this kit.
- There is no reduction in sensitivity of the COVID-19 (SARS-CoV-2) Antigen Test Midstream - Saliva against the UK variant, Brazilian variant or the South African variant. We will test the SARS-CoV-2 variants from time to time to confirm the impact of the latest variants on the test kit.

[Storage conditions & period of validity]

- Store at 4℃-30℃, and it is valid for 24 months. Do not use beyond the expiration date marked on their outer packaging.
- After the foil pouch is unsealed, the test device should be used as soon as possible and within one hour (15 ~ 30℃, Humidity ≤80%).

[Quality Control]

Program control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient volume of the specimen.

[Performance index]

- Limit of detection (LOD):** The LOD is 80 TCID₅₀/mL. A negative test result may mean that the virus concentration is lower than this value.
- High Dose Hook Effect:** When the virus concentration exceeds 2.8 x 10⁷ TCID₅₀/mL, the result may be false negative.

3. **Cross-reactivity:** There is no cross-reactivity, including human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, human coronavirus HKU1, MERS-coronavirus, SARS coronavirus, adenovirus 3, parainfluenza virus type 2, Enterovirus, respiratory syncytial virus (A), parainfluenza virus type 3, parainfluenza virus type 4a, influenza A H3N2 (Wisconsin/67/05), influenza A H1N1, influenza B (VICTORIA), Rhinovirus (HRVA30), Haemophilus influenzae, Streptococcus pneumoniae, Streptococcus pyogenes, Candida albicans, Bacillus pertussis, Mycoplasma pneumoniae, Chlamydia pneumoniae, Legionella pneumoniae, Mycobacterium tuberculosis, Pneumocystis, Pneumoniae Bacteria, human pneumonia virus (hMPV), parainfluenza virus type 1, Staphylococcus epidermidis, Streptococcus salivarius, etc.

4. **Endogenous Interference Studies:** There is no interference in studies on the following substances, including Afoxymetazoline, mucin, Budesonide, Dexamethasone, methanol, Acetylsalicylic Acid, Diphenhydramine, benzocaine, osetamivir, tobramycin, mupirocin, biotin, etc.

[Clinical Performance]

The overall study scale was 336 cases, 122 positive samples and 214 negative samples. Statistics of test results of saliva samples:

Reference RT-PCR Assay				95% Wilson Score CI		
	POS	NEG	TOTAL	PPA	NPV	Ag Test
DEEP BLUE SARS-CoV-2 Ag Test	119	213	336	97.5%	98.6%	
	122	214	336	99.5%	98.6%	
				90.9%	93.3%	
				99.2%	99.6%	
				94.9%	99.9%	
				99.2%	99.9%	
				94.1%	99.9%	
				98.6%	99.6%	

Sensitivity: 97.5% (95% CI: 90.9% - 99.2%)

Specificity: 99.5% (95% CI: 94.9% - 99.9%)

Sensitivity: Compared with the RT-PCR Assay, among people infected with SARS-CoV-2 virus, the probability of correct detection by the COVID-19 (SARS-CoV-2) Antigen Test Midstream - Saliva.

Specificity: Compared with the RT-PCR Assay, among people who have not been infected with SARS-CoV-2 virus, the probability of correct detection by the COVID-19 (SARS-CoV-2) Antigen Test Midstream - Saliva.

[References]

- Rapid SARS-CoV-2 antigen detection assay in comparison with real-time RT-PCR assay for laboratory diagnosis of COVID-19 in Thailand. *Virol J.* 2020 Nov 13;17(1):177. doi: 10.1186/s12985-020-01452-5.
- Evaluation of the Panbio COVID-19 Rapid Antigen Detection Test Device for the Screening of Patients with COVID-19. *J Clin Microbiol.* 2021 Jan 21;59(2):e02589-20. doi: 10.1128/JCM.02589-20. Prnt 2021 Jan 21.
- Detection technologies and recent developments in the diagnosis of COVID-19 infection. *Appl Microbiol Biotechnol.* 2021 Jan;105(2):441-455. doi: 10.1007/s00253-020-11061-5. Epub 2021 Jan 4.
- WHO <https://www.who.int/publications-detail/rapid-testing-for-sars-cov-2>, Visited on November 15, 2021.
- Coronaviruses. European Centre for Disease Prevention and Control. <https://www.ecdc.europa.eu/en/covid-19/latest-evidence>. Visited on November 15, 2021.

[Index/Symbols]

IVD	In vitro diagnostic medical device	⊘	Do not re-use	☀	Avoid excessive exposure to the sun
🕒	Expire date	📖	Please read the instruction for use carefully before using	📅	Date of manufacture
⚠️	Warning, please refer to the instructions in the package	🏭	Manufacturer	🚫	Don't use the product when the package is damaged
🌡️	Temperature range of product storage	📦	Batch number	⚖️	Contain sufficient quantity for $n-1$ tests
🇪🇺	European union authorization representative	☔	Keep dry	CE	CE Mark

ANHUI DEEPLUE MEDICAL TECHNOLOGY CO., LTD.
4th Floor, D-1# Zone, Pearl Industrial Park, 106 Innovation Avenue,
High-Tech Development Zone, 230088 Hefei, Anhui, China.

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Scan QR code for IFU in different languages.

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[Verweise]

- Rapid SARS-CoV-2 antigen detection assay in comparison with real-time RT-PCR assay for laboratory diagnosis of COVID-19 in Thailand. *Virol J.* 2020 Nov 13;17(1):177. doi: 10.1186/s12985-020-01452-5.
- Evaluation of the Panbio COVID-19 Rapid Antigen Detection Test Device for the Screening of Patients with COVID-19. *J Clin Microbiol.* 2021 Jan 21;59(2):e02589-20. doi: 10.1128/JCM.02589-20. Print 2021 Jan 21.
- Detection technologies and recent developments in the diagnosis of COVID-19 infection. *Appl Microbiol Biotechnol.* 2021 Jan;105(2):441-455. doi: 10.1007/s00253-020-11061-5. Epub 2021 Jan 4.
- WHO. <https://www.who.int/publications-detail/rapid-diagnostics-testing-for-sars-cov-2>. Visited on November 15, 2021.
- Coronaviruses. European Centre for Disease Prevention and Control. <https://www.ecdc.europa.eu/en/covid-19/latest-evidence>. Visited on November 15, 2021.

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ANHUI DEEPLUE MEDICAL TECHNOLOGY CO., LTD.
4th Floor, D-1# Zone, Pearl Industrial Park, 106 Innovation Avenue, High-Tech Development Zone, 230088 Hefei, Anhui, China.

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